

# The ALL-ROUND GENIUS

The three-fold solution which is as simple as it is good: Debrisoft®

Regardless of whether for deep or superficial wounds, with the Debrisoft portfolio you can rest assured of always having the right product for wound cleansing at hand.

Try Debrisoft Pad on dry, parchment-like, hyperkeratotic, seborrhoeic or oily skin.



Debrisoft® Lolly  
For deep to superficial wounds



Debrisoft® Pad\*  
13 x 20 cm / 10 x 10 cm  
For superficial wounds



**Debrisoft®**  
Debridement at its best.

\* Entire range may not yet be available.



People.Health.Care.

**Sources:**

- 1) Bahr S, Mustafi N, Hättig P, Piatkowski A, Mosti G, Reimann K, Abel M, Dini V, Restelli J, Babadagi-Hardt Z, Abbritti F, Eberlein T, Wild T, Bandl K, Schmitz M. Clinical efficacy of a new monofilament fibre-containing wound debridement product. Journal of Wound Care (2011) 20 (5): 242-248 (German translation in: Wundmanagement (2011) 5: 260-266)
- 2) Images: Barbara Pritchard, Wrexham Maelor Hospital, Wales, UK
- 3) Images: Gabriele Stern, ATZ-Stern, Bocholt, Germany

Ordering information

**Debrisoft® Pad** for debridement  
Sterile, individually sealed

Size (cm)	REF	PZN	Sales Unit (pcs)
10 x 10	34 321	131 551 58	5/50
<b>NEW</b> 13 x 20*	34 323		5/50



\* Entire range may not yet be available.

Lohmann & Rauscher GmbH & Co. KG · Westerwaldstrasse 4 · D-56579 Rengsdorf, Germany  
Lohmann & Rauscher GmbH · Johann-Schorsch-Gasse 4 · A-1141 Vienna, Austria

**AUXANO**  
MEDICAL  
AUXANO MEDICAL PTE LTD  
1 Kaki Bukit Avenue 3, #04-13 KB-1, S(416087)  
6592 5562 | www.AUXANO.com.sg


9507250 0218 e




People.Health.Care.



The Debrisoft Difference!



How Debrisoft Work?



Debrisoft.com

## Debrisoft® Pad

THE gold standard for mechanical debridement.



Debrisoft® Pad  
THE gold standard  
2 sizes, a number of possibilities

www.Lohmann-Rauscher.com


# That is **DEBRISOFT® PAD**

## The gold standard with good reason:

The patented monofilament fibres that make up the Debrisoft Pad effectively loosen wound debris, impurities and exudate from the wound, bind them to the thick fibre composite, thereby removing them effectively, yet gently from the wound. Debrisoft Pad even removes up to 99.9%\* of biofilm.

## Debridement with Debrisoft – that means:

Highly effective, yet gentle cleansing of virtually all superficial wounds.



**Debrisoft® Pad**  
For uncompromising debridement.

\* For successful biofilm management, antimicrobial treatment is required after the wound has been cleansed with Debrisoft.  
\*\* Entire range may not yet be available.

# **EFFECTIVE** as always!

Debrisoft® guarantees	Your benefits
Visible results quickly¹	You save time and money <b>Wound cleansing in 2–4 minutes</b>
Simple application¹ and immediate results	Moisten, clean and you are done <b>95% good to outstanding</b> <span>User test</span>
Virtually painless to use¹	Increased patient acceptance <b>97.4% no pain or side effects</b> <span>User test</span>



Highly efficient, effective results.

**Debrisoft® Pad**  
In a class of its own.


# **SIMPLE** for everyone!

## How to use:


**1. Moisten**  
Debrisoft Pad only cleans effectively if it is moistened beforehand! Debrisoft can be combined with everything from tap water to antimicrobial wound cleaning solutions. Whichever is selected, follow local guidelines.



**2. Clean**  
Applying gentle pressure, clean the wound using the soft fibre side. If necessary, clean the skin surrounding the wound with another moistened Debrisoft Pad.



**3. And you are done**  
Dispose of the used Debrisoft Pad. Always observe the local disposal guidelines.



**Debrisoft® Pad**  
Moisten, clean and you are done.