


Suprasorb® F Film Wound Dressing



- A** Transparent polyurethane film
- B** Bacteria-proof and waterproof
- C** Excess water vapour evaporates

Fields of application

for the treatment of non-infected superficial wounds with low levels of exudate

Indications:

- skin graft donor sites
- surgical wounds
- abrasions

Properties

- impermeable to water

Product composition

polyurethane film
polyacrylate adhesive

To note

Dressing change interval:

- To be determined by the attending physician based on the wound condition, level of exudate, and any primary dressing used, but no longer than 5-7 days.

Versions

- sterile, cut-to-size film wound dressing
- non-sterile film dressing off the roll

Suprasorb® F

Dimensions	REF	Shipping units (piece)
sterile, individually sealed		
5 cm x 7 cm	20460	10/380
5 cm x 7 cm	20461	100/1000
10 cm x 12 cm	20462	10/100
10 cm x 12 cm	20463	50/200
10 cm x 25 cm	20464	10/140
15 cm x 20 cm	20465	10/400
20 cm x 30 cm	20466	10/180
rolled, non-sterile		
10 cm x 1 m	20467	1/16
10 cm x 10 m	20468	1/16
15 cm x 10 m	20469	1/10

Not available in all countries