Rosidal® K Short Stretch Bandage

Rosidal Sys (Rosidal K + Rosidal Soft) Instructional Video





Fields of application

- for strong compression (e.g. in phlebology)
- for immobilisation of limbs
- for support and relief in traumatology and sports medicine
- for retention of dressings

Properties

- lengthwise elastic (to about 90%)
- textile elastic due to the twisting of the warp threads during weaving
- air-permeable
- textured fabric
- with selvedges
- cream-coloured

Product composition

100% cotton

To note

 can be washed up to 50 times at 95°C (commercially available mild detergent/delicates cycle)

Versions

- Rosidal K
- Rosdial Elko: 2 Rosidal K
 Bandages, 10 cm wide (for compression bandages with double bandaging technique)

Rosidal® K

Dimensions (cm x m)	REF	Shipping unit (piece)
5 m (stretched length), with bandage	clips (except for 4 cm), individually pack	ed in a folding box
4 x 5	22199	1/125
6 x 5	22200	1/108
8 x 5	22201	1/108
10 x 5	22202	1/90
5 m (stretched length), with bandage	clips, in bulk in a cardboard box	
8 x 5	22250	10/140
10 x 5	22251	10/140
12 x 5	22252	10/120
6 x 5	20650	10/90
8 x 5	20651	10/90
10 x 5	20652	10/100
12 x 5	20653	10/100
10 m (stretched length), without band	age clips, individually packed in a folding	g box
10 x 10	22206	1/60
12 x 10	22205	1/40

Elko® Rosidal K

5 m (stretched lenght), with bandage clips, pairwise packed in a folding box

Dimensions (cm x m)	REF	Shipping unit (Pack)
10 x 5	22225	1/30

Rosidal® soft Foam Bandage

Rosidal Sys (Rosidal K + Rosidal Soft) Instructional Video





Fields of application

 for padding under lymphological and phlebological compression bandages

Properties

- lengthwise and widthwise elastic
- even distribution of pressure to compensate for minor imperfections in bandaging technique (e.g. constriction)
- increases patient comfort (e.g. cushioning of the tibial crest)
- reduces slippage
- highly permeable to air and humidity
- tear-resistant and stretchable
- with cut edges
- white

Product composition

100% polyurethane

To note

- can be washed 50 times at max. 95°C
- maintains its elasticity and cushioning effect even after washing

Rosidal® soft

2.0 m unstretched length, with bandage clips, packed in pairs in a folding box

Strength (cm) pairwise in a cardboard box	Dimensions (cm x m)	REF	Shipping unit (Pack)
0.2	10 x 2	13016	1/23

Rosidal® soft

2.5 m (unstretched length), with bandage clips

Strength (cm)	Dimensions (cm x m)	REF	Shipping unit (piece)
individually packed in folding	boxes		
0.3	10.0 x 2.5	23110	1/30
0.4	10.0 x 2.5	23111	1/24
0.4	12.0 x 2.5	23112	1/16
0.4	15.0 x 2.5	23113	1/14
in bulk in a cardboard box			
0.3	10.0 x 2.5	23100	36
0.4	10.0 x 2.5	23101	24
0.4	12.0 x 2.5	23102	18
0.4	15.0 x 2.5	23103	14

Not available in all countries